

Capillary Tube *Favor COOL*



Capillary Tube

Standard according to ASTM B360

Average O.D. Tolerance: +0.03mm

Average I.D. Tolerance: +0.02mm

Condition: Soft, half hardness, hardness, ect.

O.D.		I.D.		Thickness	
(inch)	(mm)	(inch)	(mm)	(inch)	(mm)
0.072	1.83	0.026	0.660	0.023	0.584
0.072	1.83	0.028	0.711	0.022	0.558
0.081	2.06	0.031	0.787	0.025	0.635
0.081	2.06	0.033	0.838	0.024	0.609
0.087	2.21	0.039	0.991	0.024	0.609
0.093	2.36	0.042	1.07	0.0255	0.648
0.097	2.46	0.046	1.17	0.025	0.635
0.106	2.69	0.054	1.37	0.026	0.660
0.112	2.84	0.059	1.50	0.0265	0.673
0.125	3.18	0.064	1.63	0.0305	0.775
0.145	3.68	0.080	2.03	0.0325	0.826
0.188	4.78	0.130	3.30	0.029	0.737
0.200	5.08	0.145	3.68	0.0275	0.699
0.220	5.59	0.160	4.06	0.0305	0.762
0.240	6.10	0.175	4.45	0.0325	0.826

***To be adjusted as per specific requirement.**